







#### THE NEW WAY FORWARD FOR ACERS

#### WELCOME TO THE REFRESHED ACERS.

Fresh, invigorated and ready to serve our community!

#### Who is Acers?

- We believe in making your life easy, and getting things done whether its dispatching your order quickly, resolving your enquiries, actioning your product requests, or taking care to hold the right stock for you. We believe in doing it right, so its done once and it doesn't have to be done again.
- Acers is family owned, right here in New Zealand. We breathe a fresh service level and approach into dental procurement.
- The Broomhall brothers, who offer up a wide variety of skills and experiences, came together and acquired
  Acers in 2020. Backgrounds in electrical automation, stainless steel engineering and international
  consumable product development and sourcing brought about commercial skills such as quality control,
  warehousing, marketing, sales and accounting Together, this family bring together a unique range of skills
  from working at the coal-face. Experience shows that a dependable supply of quality products improves the
  end result.













#### **MISSION**

EASY, DONE.

Acers believe in making your life easy by getting things done.

#### **VALUES**

MAKE IT EASY.

GET THINGS DONE.

We're nimble, we listen and respond warmly. We do the running.

FRESH APPROACH.

We breath a fresh service level and approach into dental procurement.

FAMILY FOCUSED.

A warm family spirit that makes dentistry safe and enjoyable.

#### CONTENTS

FILLING AND RESTORATION	4
BURS	18
ENDODONTICS	25
IMPRESSION MATERIALS	28
TIPS AND NEEDLES	33
HANDPIECES / EQUIPMENT	36
XRAY	42
PREVENTATIVE AND HYGIENE	45
INFECTION CONTROL	47



# FILLING AND RESTORATIVE

#### **COMPOSITES**

PROFIL™ FLOW

#### Flowable Composite.

ProFil™ Flow is a visible - light activated, radiopaque, flowable composite. The excellent flowable consistency characteristics makes ProFil™ Flow an ideal composite for filling cavities.

#### **Features**

- · Excellent esthetic properties
- No oozing or slumping
- · Perfect for minimal invasive restorations
- · Low Viscosity
- Maximum versatility
- · High polish retention and stable color match
- · Long term fluoride release
- Low shrinkage
- · Indications
- Class III, Class V, Smaller Class IV restorations
- Base/liner in Class I and Class II restorations
- Repair of resin, porcelain and acrylic temporary materials
- Pit and Fissure sealant
- Restoration of minimally invasive cavity preparations
- Undercut blockout.



Sensitivity to ambient light: Above 100s.

Weight of total inorganic filler: 78%

· Volume of total inorganic filler: 59%

Depth of cure: 1.5 mm

Compressive strength: 230 MPa

Flexural Strength: 120MPa

• Water Absorption: 40 µg / Щ

Solubility: 7.5 µg / Щ.



8gr, Advanced Resin Technolo 63A01 Profil Flow A1 63A02 Profil Flow A2 OUR PRICE

AL OR





#### PROFIL MICRO HYBRID COMPOSITE

Universal Restorative Composite Resin, for use in both Anteriors and Posteriors. ProFil™, is a micro hybrid composite material available in Syringe or Capsules.



#### **Features**

- · High polishability
- Non sticky formula
- · Optimal balance of physical properties
- · Light cured, radiopaque and color stable
- · Low polymerization shrinkage
- Fast placement technique, easy to apply and sculpt
- · High strength and wear resistance
- · Available in 14 shade.

#### **Physical Properties**

ProFil™, a micro hybrid composite has outstanding physical properties that provide excellent wear resistance.

- Sensitivity to ambient light Above 100 sec.
- · Weight of total inorganic filler 78%
- · Volume of total inorganic filler 59%
- Depth of cure > 1.5 mm
- · Compressive strength > 230 MPa
- Flexural Strength > 120MPa
- Water Absorption < 40 µg / Щ
- Solubility < 7.5 µg / Щ.

#### 20x.25g Compules



4g Syringe



20 x 0.25gr, 33% Extra Free Material, P, Micro Hybrid, 18CA01 Profil Capsules A1 18CA02 Profil Capsules A2 18CA03 Profil Capsules A3

18CB03 Profil Capsules B3

OUR \$69.00

. 184A01 Profil Syringe A1 184A02 Profil Syringe A2

#### PRIME-DENT® FLOWABLE COMPOSITE

- Highly polishable, 67% filled esthetic composite used for all types of cavity preparations.
- · Kit includes 4 x 2g syringes in shades A2, A3, A3, 5, B2,













#### PRIME-DENT® NANO-HYBRID COMPOSITE



A resin composite with unique nanofiller technology. Formulated with nanocluster filler particles, it combines the strength of a Hybrid with the polish of a Microfil.

- · Easy to control during placement.
- · Not sensitive to ambient light.
- · Adapts well to cavity preparations.
- · Does not stick to the placement instruments.
- Can be contoured to the desired shape before curing.
- Highly rated for blending with tooth structure, radiopacity and ease of finishing.
- Finished and polished to a high luster.



Kit: 001-610 OUR PRICE
Prime-Dent Nano Hybrid
Composite Kit - Surinae

\$124.00



#### **EZETCH**

#### 37% Phosphoric Acid Etchant.

EZ Etch contains phosphoric acid for etching tooth enamel and dentin when preparing surfaces for application of composites.

- High quality, non drying, medium viscosity gel
- · Stays put, does not run
- Large & comfortable to use syringes
- · Visible color
- · Non drying
- · Washes off easily

- · Easily placed
- · One lump wash off
- 5x3ml syringes + 10 tips
- Tip can be angled to facilitate easy and accurate placement of the gel even in disal and lingual composite preparations











#### ECO-STM

#### Opaque Resin-based Pit and Fissure Sealant.

- Eco-S<sup>™</sup> is a single part, light curing, pit and fissure sealant.
- It is Based on methacrylate resin chemistry. It can be used as a pit and fissure sealant for the prevention of caries.
- Curing Time 20 seconds.
- Pack of 4 x 1.2ml syringes with 20 disposable prebent tips.



#### Features and benefits

Eco- $S^{\mathsf{TM}}$ , pit and fissure sealant, isalight - cured dental material designed for sealing the enamel pit and fissure of teeth for the prevention of caries. By adding specific nano-fillers to Eco- $S^{\mathsf{TM}}$ , it has good properties like shear bond strength, micro-leakage and wear resistance.

 $Eco-S^{\mathsf{TM}}$  has an appropriate viscosity, to enable the sealant to penetrate the depths of pits. So it provides good efficacy as well as workability.

#### **Indications**

Sealing the enamel pits and fissures of teeth for prevention of caries.

#### Package

 $4 \times 1.2 \text{ml}$  syringes with 20 disposable prebent tips.









#### **PRIMA 2000**

#### Single Step Bonding Agent

- · Primer with adhesive
- Use with Autocure
   Activator to enhance
   bond in light inaccessible
   procedures
- 10ml bottle





#### **PROLINK**

#### **Universal Dental Adhesive**

- · Primer and sealer in one bottle
- · Minimal treatment time
- Perfect marginal integrity
- · Reduced polymerisation shrinkage
- Increased adhesion & durability
- 5ml bottle





#### **CEMENT**

#### HIGH Q BOND SE TRANSLUCENT

#### Self-Etch/Self-Adhesive Dual Cured. Permament Adhesive Resin Cement.

- Superior retention and total margin integrity
- Does not require etching, priming or bonding before cementing the permanent restoration.
- Self-cured cement which you can lightcure for an immediate set and extra strength at the margins.
- Radiopaque Easily seen in radiograph.
- Fluoride Known to reduce the possibility of secondary caries and considered to be cariostatic.
- Auto mix Syringe Saves application time, guarantees consistent mix.
- Virtually neither taste nor odour Increased patient comfort.











#### **AMALGAM**

#### PRO COPPER

#### Amalgam 40% Silver

- · Zinc free
- Positive dimensional change assisting marginal sealing
- · Excellent handling
- · Uniform consistency
- 50 capsules/box.



OUR \$119.00

OUR \$93.00



#### **PRO SILVER**

#### Amalgam 67.9% Silver

- · Zinc free
- High content of silver and copper evenly dispersed
- 50 capsules/box.



OUR \$199.00




#### **MATRICES**

#### SADDLE CONTOURED METAL MATRICES UNIVERSAL KIT

#### Shape with enlarged undergingiva ledge.

Edges of saddle matrices have cylindershape tubes where tines or a spring clip are inserted upon matrix installation.

- · Easy-to-use matrix system
- Optimal approximal contact without difficulty
- · Three lengths
- · Stainless steel
- 0.050mm thick.











#### METAL CONTOURED MATRICES

- Anatomically pre-shaped
- Shapes for molars and premolars

Provide perfect approximal form and contact

- · Stainless steel
- 0.050mm.



without ledge





For Premolars bilateral





For Molars without ledge





For Molars bilateral



Matrices 1.501 and, 1.511 are mainly to be used with Tofflemire retainer.

#### TRANSPARENT CONTOURED MATRICES KIT

- Anatomically pre-shaped.
- Four shapes of premolar matrices and four shapes of molar matrices.
- Provide perfect approximal form and contact.





#### For premolars:

shape 1 (one central ledge) 15 pcs.



shape 2 (two central ledges) 15 pcs.





shape 3 (right ledge) 15 pcs.



shape 4 (left ledge) 15 pcs.

#### For molars:



shape 1 (one central ledge) 15 pcs.



shape 2 (two central ledges) 15 pcs.



shape 3 (right ledge) 15 pcs.



shape 4 (left ledge) 15 pcs.

#### TRANSPARENT SECTIONAL CONTOURED MATRICES















#### **ANTERIOR MATRICES & STRIPS KIT**



Proximal anterior strips can be used upon restoration of proximal surfaces of anteriors.



Proximal anterior matrices 1.387A – 6 pcs.




Proximal anterior matrices 1.387 – 6 pcs.







Twin anterior matrices provide







# SECTIONAL CONTOURED METAL MATRICES

Matrices of 0.05mm thickness are intended for restoration of large decays in subgingival area. Such matrices retain the shape well.

- 0.050mm, Hard
- Anatomically pre-shaped
- Provide tight anatomical contact
- · Stainless steel
- Hard matrices are easy to install
- 50pcs





Medium



Large



Large with Ledge

#### METAL MATRIX BANDS

#### For Premolars

Finest Tofflemire bands available.

Consistently achieves better contours, tighter contacts, and healthy periodontium in posterior composite or amalgam restorations with easy to apply matrix bands.

#### For Molars

#### Without Ledge

- Finest Tofflemire bands available.
- Consistently achieves better contours, tighter contacts, and healthy periodontium in posterior composite or amalgam restorations with easy to apply matrix bands.
- 0.050mm
- 12pcs













#### FIXING WOODEN WEDGES



Superthin Short





Thin Short





**Medium Short** 





Thick Long



#### TRANSPARENT PLASTIC **WEDGES KIT**

- For Interdental Restorations 80pc
- Thin (1.5 x 1.5 mm) 50pcs
- Medium (1.5 x 2.0 mm) -30pcs.





#### **ELASTIC FIXING WEDGES KIT**

- Elastic Fixing Wedges Kit 40 pcs.
- 2.0 mm thick x 20 pcs.
- 2.5 mm thick x 20 pcs.






#### **STANDARD** RING

- · Tines height 5 mm.
- · High-Quality Stainless Steel Clamp with Flat Tines to be used with Forceps.

#### M D RING



· High-Quality Stainless Steel Clamp with Flat Tines to be used with Forceps.











#### **DELTA RING**



- · Delta Ring's double tines are suitable for the installation of fixing wedges and add-on wedges.
- · High-Quality Stainless Steel Clamp with Flat Tines to be used with Forceps.











#### **FORCEPS**

#### For installation of sectional matrices.





#### METAL STRIPROLL

- 6mm x 3m
- 0.050mm thick.

Code: 1.392
OUR \$30.00



#### TRANSPARENT STRIPS

- 10 mm wide x 95 mm long
- 0.05mm thick 100pcs
- Transparent anterior strips
- No hindrance to light-curing.



Code: 1.041

OUR \$11.00

# TRANSPARENT STRIPROLL WITH CUTTER

- 10m long x 10mm wide
- Especially indicated for light-cured composite materials.
- The matrix strip is highly resistant and does not tear under tension – In square box with cutter.



# • 10m long x 10mm wide

TRANSPARENT STRIPROLL

- Tom long x Tomm wide
   Continuous roll matrix
   band made of transparent
   polyester, especially
   indicated for light-cured
   composite materials.
- The matrix strip is highly resistant and does not tear under tension.



Code:1.245

OUR \$24.00

Code:1.241

OUR \$13.00







#### TRANSPARENT CROWNS UNIVERSAL KIT

Single Transparent Crowns available in Anterior, Posterior, or Deciduous kits. Individual shapes for incisors, canines, premolars and molars.



#### Anterior Transparent Crowns Kit

- Anterior Transparent Crowns Kit 64 pcs.
- 44 pcs for incisors (upper & lower jaw, small, medium, large).
- 20 pcs for canines (upper & lower jaw, small, medium, large).







#### Posterior Transparent Crowns Kit

- Posterior Transparent Crowns Kit 64 pcs.
- 32 pcs for premolars (upper & lower jaw, small, medium, large).
- 32 pcs for molars (upper & lower jaw, small, medium, large).







#### Deciduous Transparent Crowns Kit

- For premolars and molars (upper and lower jaw, small, medium, large) - 32 pcs.
- For incisors and cannines (upper and lower jaw, small, medium, large) - 32 pcs.
- Crowns are available in three sizes: small (S), medium (M), large (L).











#### **POLISHING**

#### PROFIL FINISHING KIT

#### **Autoclavable Finishing Burs**



- The set is designed to make the tough task of composite finishing significantly simpler, shorter, and more efficient while giving excellent clinical results.
- Polish the composite without the need to use disc, paste, or polisher, with the benefit of a reusable autoclavable solution that saves cost over time.
- No need to use polishing discs/rubber polishers
- Faster work 10-20 seconds instead of 3-5 minutes
- Much more cost effective save up to 75%
- · No bleeding doesn't cut the gingival
- · Quick and easy learning curve. Very durable.

Code: 199PFK
OUR \$90.00

#### WHITE STONE POLISHING KIT

- FG White Stone Point composite finishing Kit/12.
- Made of micro-grained aluminium oxide grains.
- Excellent for fine finishing of composite restorations, compomers, enamel and porcelain. Contains Flame (F4) x 6 and Needle (F5) x 6.

Code: WSPFK
OUR \$36.00



#### **MANDRELS**

#### Polishing discs



Stem Polishina Discs: 1.121



Stem Polishing Discs: 1.121
OUR \$19.00

#### **TEMPORARY**

#### **PROFIL TEMP**

- · Temporary Filling Material
- Similar to Cavit
- Rapid initial setting
- Good adhesion to dentin
- · Available in 28g jar



Code: 19T028

OUR \$4.7.00

#### QTEMP



- Temporary Cement. Not irritating to soft tissue will not inhibit permanent cementation
- Two-Stage Curing Initial gel-set stabilizes temporary, allowing easy removal of excess.
- Automix syringe. 5ml syringe.









#### **POLISHING DISCS**

- · Easy-to-use finishing and polishing system allows you to create high-luster restorations that offer life like esthetics.
- Discs are color coded to allow for easy identification of grit sequence
- · Flexible disc technology makes changing discs easier to save you time
- · Discs conform to the contours of the teeth for easier finishing and polishing
- Extra thin discs allow for easier access to interproximal surfaces while offering a slightly stiffer disc.



#### Stem Polishing Discs





OUR \$34.00

#### Polishing Discs 12mm & 16mm diameter



Polishing Discs Fine For Finishing Fine (White) 12mm Code: 1.733-12 Fine (White) 16mm Code: 1.733-16

Polishing Discs Medium For Contouring Medium (Green) 12mm Code: 1.732-12 Medium (Green) 16mm Code: 1.732-16

Polishing Discs Coarse For Gross Reduction Coarse (Blue) 12mm Code: 1.731-12 Coarse (Blue) 16mm Code: 1.731-16

OUR \$36.00

#### **POLISHING STRIPS**

Fine For Finishing/Polishing (White) - Code: 1.052

Medium For Contouring (Green) - Code: 1.051

Coarse For Gross Reduction/Contouring (Blue) - Code: 1.050



OUR \$12.00





#### **ACCESSORIES**

#### **COTTON ROLL HOLDER**

· Cotton roll holder for fixation of cotton rolls during dental procedures.



OUR \$19.90

#### MIXING PADS

- Non-slip
- 50 sheets
- Three sizes: S: 5cmx4cm L: 7.5cmx5.5cm XL: 13.5cmx9.5cm.



**\$7.00** 





#### MICRO APPLICATORS

- Tube package contains 100pcs × 4tubes/1 box = 400 applicators.
- · Multiple applicator colors are available to identify different materials during a procedure
- Non-absorbent, non-linting burrs hold a drop of solution in suspension without dripping or spilling
- · Micro applicators can be used for bonding agents, cavity linings, sealants, etchants, dycal applications, uride varnishing, hemostatic solutions, disclosing solutions.



#### **BRUSH TIP** AND HOLDER

· White color only.













#### **FG DIAMOND BURS**




Diamonds are engineered for multi-use and made from 100% natural diamonds.

Shanks are made from stainless steel. Shank hardness 40HRC. High precision computerized manufacturing process.

All diamonds guarantee 100% satisfaction.

\$42.00


#### **ROUND HEAD**








#### **INVERTED CONE**



#### INVERTED CONE W/ COLLAR CYLINDERS ROUND END







#### **PEARS**











#### **CYLINDERS**



#### **ROUND END TAPERS**




#### POINTED CONE / NEEDLE




PACKET OF 10 \$**42.00** 







#### FLAT END TAPERS



#### **OGIVAL END TAPERS**

Code	DBSF878K/12	DBSF878K/16	DBF878K/12	DB878K/10	DB878K/12	DB878K/14	DB878K/1
Grit	Super fine	Super fine	Fine	Medium	Medium	Medium	Medium
Head size	1.2	1.6	1.2	1.0	1.2	1.4	1.6
Head Length	8.0	8.0	8.0	8.0	8.0	8.0	8.0
Total length	21.0	21.0	21.0	21.0	21.0	21.0	21.0







#### **TURBO DIAMOND BURS**

### The design of the Turbo Diamond Bur provides a special cooling effect to prevent heat increase.

These improved second-generation Burs will work faster for you at a cooler temperature.

- The Turbo Bur design rapidly eliminates excess heat to allow a cooler operating temperature reducing the level of tissue trauma
- · Bur design enables water through the spiral channels directly to the contact point
- Self-cleaning for faster removal of debris. Reduces the danger of pulp lesions
- Multi-layer cutting edge with spiral channels coated in natural diamond particles to support faster, super high cutting efficiency.

Code	TDB850/014	TDB856/014	TDBS850/018	TDBS856/018
Grit	Coarse	Coarse	Coarse	Coarse
Head size	1.4	1.6	1.8	1.8
<b>Head Length</b>	10.0	10.0	10.0	8.0
Total length	23.0	23.0	21.0	21.0
		'	Short	Short

\$75.00

#### **TURBO CROWN REMOVAL BURS**

#### Revolutionary Crown cutting solution.

Hybrid Twister water cooling crown cutter, the Turbo Crown bur is designed to operate faster at a cooler temperature by delivering water to the contact point. Designed for Zirconia crown, and metal crown or bridge removal. No breakage, non slip, no vibration.

- The Turbo Bur design enables water through the spiral channels/grooves directly to the contact point
- Rapidly eliminates excessive heat to allow a cooler operating temperature reducing the level of tissue trauma. Multi-layer cutting edge
- The spiral channels/grooves are all coated in natural diamond particles to support super high cutting efficiency and fast removal of debris.

-		
1	Code	CT1DE
ALC:	Grit	Coars
19	Head size	1.8
442	Head Length	6.0
185	Total length	21.0
1980		

\$**87.00** 









#### **FG SURGICAL BURS**

#### DIAMOND ORAL SURGICAL NEEDLE

	A
Code	DBOS52L
Head size	1.2
Head Length	10.0
Total length	32.0


\$59.00

#### SUPER MINI DIAMOND ORAL SURGICAL ROUND END TAPER



PACKET OF 10 \$**59.00** 

#### CARBIDE ORAL SURGICAL ROUND HEAD



PACKET OF 10 \$**59.00** 

#### **CARBIDE SURGICAL 151**

• Cutting wisdom teeth and accessing deep tooth root.

Code	CB151
Head size	1.4
Head Length	11.0
Total length	25.0

PACKET OF 3 \$49.00







#### **FG SURGICAL BURS**

#### **CARBIDE SURGICAL 2158**



- Cutting molars or wisdom teeth, especially when the patients mouth is small or narrow.
- CB2158 can be used to cut a groove in the molar or wisdom teeth, and CB151 used to cut the tooth if the tooth is deep.

Code	CB2158
Head size	1.2
Head Length	4.0
Total length	19.0



#### DIAMOND ORAL SURGICAL BUR ROUND



 Cutting wisdom teeth, such as: distoangular impaction, vertical impaction, mesioangular impaction and horizontal impaction. Integrated one piece and will not break during operation. Can also solve depth wisdom teeth cutting issue.

Code	DBOSD14	DBOSD16
Head size	1.4	1.6
Head Length	1.4	1.6
Total length	25.0	25.0

PACKET OF 10 \$63.00

#### FG DIAMOND ENDODONTIC BUR

• Minimally invasive instrument for root canal treatment. Easily access pulp chamber. Excellent treatment sight.

Code DBEN1		
Head size	1.5	
<b>Head Length</b>	10.0	
Total length	22.0	



# ENDODONTICS

#### WELL-ROOT™ PT

## Premixed root canal filling material



#### Features and benefits

- · Premixed putty without mixing process
- · Paste-only filling
- Proper setting time for applications
- · Superior compressive strength
- · Antibacterial effects by alkaline pH
- Biocompatible elements
- · Radiopaque.
- Calcium aluminosilicate paste (putty type)
- 0.25g x 10 capsules

#### **Indications**

- Pulp capping
- Repair of root perforation
- Repair of root resorption
- Root-end filling
- Apexification.

Code: 2070
OUR \$174.00



#### Properties measured according to ISO 6876:2012

Flow (mm)	Film thickness (un)	Setting t	ime (min.)	Dardin and the Assess Al
Flow (mm)	Film thickness (an)	Initial	Final	Radiopacity(mmAl
9.3±0.2	38.8±4	5±2min	45±5min	5.4







#### WELL-ROOT™ ST

#### Premixed root canal filling material

"Good product, value for money and easy to use."

Stephen Shepherd BDS East Bay Dental Centre



- Premixed, injectable, bioactive, calcium aluminosilicate compound paste.
- For permanent obturation of the Root Canal
- Conveniently premixed ready-to use
- · Requires the presence of water to set and harden
- · Does not shrink during setting
- 1 x 2g syringe + disposable tips + holder.



Code: 2077

OUR \$129.00

#### **GUTTA PERCHA**

- Round Full Kits #15-50
- 240 points





#### PAPER POINTS

- Round Kit #15-50
- 400 points

Code: 2062
OUR \$19.00







#### **WELL-PASTE**



#### **WELL-PEX**







#### **Temporary Root** Canal Filling 2g

- · Filling material for infected root canal
- Contains calcium hydroxide and barium sulfate
- · Highly flowable premixed paste.

#### **EDTA Root Canal Preparation Cream**

- · Easy and fast removal of debris
- Facile cleansing and shaping of the root canal
- Lubricating effect simplifies its injection into canals
- 3x6g syringes, disposable tips.

#### Oil Based Root **Canal Filling**

- · Filling material for infected root canal
- · Contains calcium hydroxide and iodoform
- · Good flowability and antibacterial characteristics
- Contains 2g syringe, disposable tips and extension holder.

#### **ENDO CHUCK INSERT**

- The Bonart endo chuck insert is designed to support endodontic procedures.
  - The end of the tip has a slot to hold a U-file.
- · U-file not included.
- For use with 25KHz magneto ultrasonic scalers only.
- CE and ISO certified 510K registered
- Tip: Stainless Steel
- · Stem: Stacked nickel strip.



OUR \$150.00







# IMPRESSION MATERIALS

#### **VPS**

#### 3 HOLE DISPENSING GUN HOLDER

- Easy to organise guns/tips
- Transparent.









#### HYDROPHILIC VINYL POLY SILOXANE



- · Highly hydrophilic
- Exceptional accuracy in recording soft and hard tissues
- · High tear strength
- · Excellent dimensional stability
- · Optimal elastic recovery
- · Snap-set technology
- Cost-effective (4 cartridges in each box).



#### **AQR8 LIGHT**

#### Light body vps



 Regular Set Setting time:
 3.5 minutes



- Fast SetSetting time:2 minutes
- Maximum precision
- High-tear strength
- Excellent performance in bleeding spots.

Regular: 3003 | Fast: 3008

OUR \$95.00

#### **AQR8 MONO**

#### Monophase Medium vps



- Working time: 2 minutes.
- Setting time:3.5 minutes.
- Maximum precision
- High-tear strength
- Excellent elastic recovery.

Medium: 3004

OUR \$95.00

#### AQR8 HEAVY

#### Heavy body vps



Regular SetSetting time:3.5 minutes



Fast Set
Setting time:
2 minutes

#### AQR8 BITE REGISTRATION FAST SET



- Easy handling
- Scannable
- · Rigid and highly stable
- Very accurate
- 4 x 50ml cartridges

Code: 3005
OUR \$95.00







#### **DENTIANN™ HYGEDENT®**

#### Alginate Impression Material, Regular Set.

- · DentiAnn Alginate Regular Set Mint 454g bag
- Great definition and detail
- · High Elasticity
- · Nice mix and texture
- · Improved tear resistance
- · Mint flavour
- · Dust free
- · Mixing time: 30 secs.
- · Working time: 1 min 30 secs.
- · Time in mouth: 1 min 30 secs.
- · Setting time: 3 mins.
- · 22 complete impressions.





#### **EAZIPAN LC**

Eazipan LC is a light-cured transparent thixotropic impression tray material in form of coloured plates for precise work for crown, bridge and total prosthetic technique as well as for fabrication of base for bite patterns.

#### **Advantages**

- · Easy handling
- No offensive smell
- · No toxic fumes and waste
- · Minimal shrinkage
- Each box contains 50 sheets.

#### **Indications**

- · Individual tray for use in crown, bridge and total protestic technique
- · Base plates for bite patterns and odontograms.







#### **IMPRESSION TRAYS**

- · Disposable and sanitary
- Cost effective
- · Shape adjustable with any heating system
- · No tray cleaner or adhesive needed
- Semi-transparent tray for easy viewing and checking
- · Echelon-style holes improving hold on impression material
- · With denture bur alter the tray to use it for implant.
- 10 different sizes: Large: Upper, Lower Medium: Upper, Lower Small: Upper, Lower Partial: Left, Right Anterior: Upper, Lower



#### **▶ UPPER TRAY**

	w	н
SIZE(mm) DT-UL	78	61
DT-UM	78	61
DT-US	78	61

#### **▶ LOWER TRAY**

	W	Н
SIZE(mm) DT-LL	84	61
DT-LM	77	54
DT-LS	72	51

#### ► PARTIAL TRAY (LEFT/RIGHT)

	W	H
SIZE(mm) PB	51	53
PA	51	53
P1	70	38
P2	66	38

#### OUR \$48.00



**Upper Large** 40/baa Code: SGTRAYLU



**Upper Med** 48/bag Code: SGTRAYMU



**Upper Small** . 48/baa Code: SGTRAYSU



Quadrant 48/bag Code: SGTRAYPB



Quadrant 48/baa Code: SGTRAYPA



**Lower Large** 40/bag Code: SGTRAYLL



Lower Med 48/bag Code: SGTRAYML



Lower Small 48/bag Code: SGTRAYSL



**Upper Anterior** 48/bag Code: SGTRAYPAU



Lower Anterior 48/bag Code: SGTRAYPAL







#### MATRIX BUTTONS



- Make custom impressions for temporary crowns
- · Soften in hot water then form around teeth
- In and out of the patients mouth in 30 seconds
- Additional uses: insulate sensitive teeth while operating in another area.
- Implant shims to guide placement.

Code: 2073
OUR \$79.00



- · Short setting time in the mouth
- Non-greasy, smooth, easy to mix
- No air bubbles or air traps
- · Very good mechanical properties
- Excellent reproduction of details
- · Hydrophilic.

Code: 3007
OUR \$94.00

#### **BITE TRAYS**

#### Disposable • Perfect impression result • No contraction • Ultra thin contact



Posterior 50/box Code: SGMULTIP

OUR \$26.00



Anterior 30/box Code: SGMULTIA

OUR \$26.00



Quadrant 30/box Code: SGMULTIQ

OUR \$26.00

# NEEDLE AND TIPS



- Special design for effectiveness of sterilization
- Electric discharge of all products on hub hole
- Tri- bevel point for easy tissue penetration without coring
- · The round cutting edge protecting the gum
- The most strict grinding angle, ratio of tri-bevel points reduce pain
- · The highest quality of saline coating substantially reduces pain and smooth repeated injection
- Special heat treatment of pipes makes super-elastic
- Bevel indicating mark.

\$19.00







#### MAIN MIXING TIPS



- Yellow Mixing Tip S130: Light body impression material.
- · Neo Green Mixing Tip - S131 Heavy body impression material, Bite registration.
- · Neo Pink Mixing Tip -S132: Use for regular body impression material.
- · Neo Blue Mixing Tip -S133: Core build-up material, Light body impression material.
- · Neo Core Mixing Tip - S123: for build-up material.
- · Short Resin Mixing Tip - S124: Core build-up material.
- · Long Resin Mixing Tip - S125: Core build-up material.
- Blue-Orange Mixing Tip - S176. Irrigation and injection into root canals.
- · Blue-White Mixing Tip - S177. Irrigation and injection into root canals.


OUR \$26.00

#### **ACCESSORIES**

- White intraoral - S116: Regular body impression material.
- Core Intra Oral Tip (S117): for Neo Core Mixing Tip.
- Core Intra Oral Tip - L (S118): for Neo Core Mix Tip-L.
- · Capillary Endo -S121: Irrigation and injection into root canals.
- Bite Intra Oral Tip (S122): for Neo Green Mixing Tip. Heavy body impression material, Bite registration.
- Yellow Intra Oral Tip (S115): for Neo Yellow Mixing Tip and Neo Blue Mixing Tip.



#### AUTO **MIXING TIPS**



Yellow Neo Auto Mixing Tip (S135): for Automatic mixing unit (Heraeus Kerr, Zemack Dentron Omicron Draue).

OUR \$34.00



#### INTRAORAL **NEEDLE TIPS**



- Brown for core mixing tip
- Yellow for yellow mixing tip
- Easily bent
- 50 tips/pack





# DISPOSABLE SYRINGES







#### **DENTAL INTRAORAL TIPS**

 For use with Amalgam and GIC capsules for more accurate/delicate placement of product.





#### MIRACLE TIP

- · Low Viscosity Tip
- Easily fill the compule with your restorative material such as GIC mixes, put the tiny plug into the compule and load it into your composite gun
- Easy to inject into post preps
- Compatible with all C-R syringes and Resin Guns.







#### MIRACLE TIP WITH NEEDLE

- · Easily bendable stainless-steel needle tip
- Easily fill the compule with your restorative material such as GIC mixes, put the tiny plug into the compule and load it into your composite gun.
- · Perfect for flowable composites, core material, cements, and glass ionomers.
- Easily inject anything into the tooth or restorations.





#### **PREBENT TIPS**



- Black 20 Gauge Tip: for Flowable Composite.
- Black 25 Gauge Tip: for Etch or Flowable Composite.
- Blue 25 Gauge Tip: for Etch.

Blue 25G: 4052 | Black 20G: 4050 | Black 25G: 4051

OUR \$20.00



# HANDPIECES/ EQUIPMENT

#### **HIGH-SPEED HANDPIECES**

**MAXICURE TIP** 





#### **MAXICUT**









- For CA bur (ø2.35).
   Ball bearing cartridge
- Max speed 30,000 min.



#### MAXI-E-MOTOR LED

#### Brushless Electrical Micro Motor

- Speed 2,000~40,000rpm. Power: 40w
- · Compatible with all E type handpieces
- Includes: motor, cable, tube, circuit board and foot control plug.











#### **MAXISTRAIGHT**

#### Fibre Optic Straight Hand piece with Internal Spray

- For HP burs (ø2.35)
- E-type connector
- · Internal water spray
- · Speed: 22,000-27,000rpm.







#### **MAXITORQUE**

#### 45° Surgical Fibre Optic



- Kavo, NSK and W&H connection options available
- · New design of in-air system better airflow, gruntier power, and quieter
- · Triple water spray for full cooling
- · Bright light excellent visibility
- GC German Ceramic bearings
- Head size 11.0mm wide x 13.3mm high
- · Fibre-optic light
- Colour temperature (K) 5,900
- Brightness (LUX) 25,500
- Mpa 0.25-0.3
- 22 watt
- Speed up to 300,000 rpm
- Dynamic Balancing Correction
- Steam sterilisable up to 135°C
- Kavo connection, NSK connection, or W&H connection options available.











#### **MAXITORQUE**

#### High Speed Kavo Fibre Optic Air Turbine

- · Kavo connection
- · Triple water spray
- Steam sterilisable to 135°C
- Speed 280,000rpm
- 25 watt MPa 0.25-0.3
- · Dynamic balancing correction.









#### **MAXISPEED**

#### High Speed NSK Mach connection Fibre Optic Air Turbine

- · Standard NSK Mach connection
- GC German Ceramic bearings
- Super Torque mini head size 11.0mm wide x 13.3mm high
- · Fibre-optic light
- Colour temperature (K) 5,500
- Brightness (LUX) 25,500
- Mpa 0.25-0.3
- · Triple Port water spray
- 22 watt
- Speed up to 300,000 rpm
- · Dynamic Balancing Correction
- Steam sterilisable up to 135°C.













#### **MAXIPACE**

#### High Speed W&H connection Air Turbine with LED

- · W&H connection
- GC German Ceramic bearings
- Super Torque head size 12.7mm wide x 13.1mm high
- · LED light
- Colour temperature (K) 5,500
- Brightness (LUX) 25,500
- Mpa 0.25-0.3
- · Triple Port water spray
- 25 watt
- Speed up to 280,000 rpm
- · Dynamic Balancing Correction
- Steam sterilisable up to 135°C









#### C-MORE HD MIRROR


#### High definition, high reflectance mouth mirror

- · More reflective light seem more clinically!
- Superior reflection, sharpness, and contrast without distortion
- Find root canals, margins, and fine details more easily.









C-More Mirror







#### MIRROR SUCTION SYSTEM

- Dentiann 5-Mirror Suction System
- Four-handed Dentistry? No! Only two-handed. Mirror and Suction at the same time.
- · Easy to operate without assistant.
- Kit includes 5 x Anti-Fog Mirrors + 2 x Silicone Hoses + 3 x Adapters +1 Mirror Cleaning Brush.
- Mirrors: Front mirroring. Anti-Fogging function.
- Suction: 8 holes 4 holes on the side of mirror and 4 holes on the bottom of mirror.
- · Excellent suction function.
- · Adapters: Install to any kind of unit with 11mm, 16mm, and 18mm adapters.
- Autoclavable: mirror heads, handle, and silicone hoses are all autoclavable.








- Front mirroring
- Anti-fogging function REF S901-1


6 x Mirror Suction Heads



#### **High compatibility**

Install to any kind of Unit Chair with 11mm, 16mm, 18mm adapters








#### **8 Suction Holes**

- 4 holes on the side of the mirror and 4 holes
- on the bottom of the mirror - Excellent suction function
- 8 suction holes do not suck tongue or any parts inside of the mouth.



#### Easy to clean

- Cleaning the inside of the mirror with cleaning brush
- Mirror head, handle and sili cone hose are all autoclavable



#### Simple Installation

Easily Connect to your current HVE valve









#### **CUREBOX**

#### **Light Curing Unit**

- · Protected from overcurrent and overheating
- Stable power supply
- Faster curing using an LED lamp 4 times stronger than before
- · Selectable wavelength
- · Adjustable brightness.

OUR \$449,00



#### **DENJOY LED CURING LIGHT 7W**



- Wave length 400-515nm blue light spectrum.
- High output, curing 3.5mm resin within 10 seconds.
- · 2000-2400 mw/cm2.
- At full charge can be used 200 times with a 10s curing time.
- · LCD Display.
- Curing Time 5/10/15/20 seconds. Beeps every 5sec.
- · Three curing programs Full/Ramp/Ramp to half
- Available in White/Orange colour.



#### **AMALGAMATOR**



- The universal high-speed amalgamator is controlled by a microprocessor, which ensures precise, controllable and consistent mixing.
- · With a safety cover, the amalgamator will deactivate automatically when you open the front cover to access the capsule.
- · Microproccessor ensures precise, controlled consistent mix
- Adjust speed from 2800-5000rpm
- Digital LCD displays time and speed
- Quiet, stable operation
- Motor stops when mask is open
- Available in orange.














- · Use for one tooth or an entire quadrant
- Good on wet or dry teeth
- Good on metal and non-metal surfaces
- 13 microns.



#### SENSOR WIRE CABLE SAVER

- · Carts, feet, tangles, can damage your digital sensor wire.
- New Cable Saver wraps your sensor wire in a cusion of protection. It's color coded so you can easily avoid accidents too!
- One size fits all!
- Available in 8inch lengths.
- Simply cut to size!





#### COMFEE'S ECONO


#### Sensor / PSP Protective cover

- · Protected against tears
- Super-soft vinyl material
- Rounded corners
- Unique tear-away seams.



#### **SNUGGIES** SENSOR SLEEVES

Works with sensors of any size, shape, and thickness one-size-fits-all.



Simply peel, stick, and go.

#### **BITE BLOCK BAGGIES**



- Bite Block Baggies protect your equipment as well as your patients
- One bag covers film, PSP, as well as the sensor and alignment arm!
- Handy selfdispensing box.







#### PSP PROTECTIVE COVER / BARRIER ENVELOPES

#### **PSP Protective Cover**

#### ECONO

· With thumb notch for unloading

 Black cover stops light exposure.

Size 0 - 80100 Size 1 - 80110 Size 2 - 80120 Size 3 - 80130

\$89.00



#### **DELUXE**

• Easy-Tear for quick and easy plate removal

 Large Tabs on protective covers facilitates easy sealing of the envelope when wearing gloves.

Size 0 - 80101 Size 1 - 80111 Size 2 - 80121 Size 3 - 80131

\$99.00



#### OPT

- Pre-loaded protective covers with cardboard inserts for plates with magnets, makes loading effortless while eliminating the possibility of accidentally loading the plate backwards
- Large Tabs on protective covers facilitates easy sealing of the envelope when wearing gloves
- Handy storage box keeps envelopes neatly stacked for faster loading and less fumbling.





\$**120**.00













#### STICKEM'S™



- Packed in automatic dispenser box
- · Simply pull, fold and stick.

OUR \$50.00



#### STICKY **BITES**



- X-Ray Positioning Foam Fast and easy positioning.
- Soft polyethylene foam
- For Anteriors, Posteriors or BWs. 100% latex free.

OUR \$50.00



#### DIGI LOOPS



- · X-Ray Sensor Heavy duty paperboard
- Colour-coded for sizes.



OUR \$31.00

#### SUPA'S

- · Foam Film/ PSP Holder
- · Firm holding at 90 degrees
- · For perfect alignment use with RAPiD.







#### **BITF-WING** LOOPS

- Paper Loop Film/ PSP/Sensor Holder
- Made of heavy duty paper board
- Great for taking bite wings.



#### **CLEAN UP'S**

- **PSP Cleaning Wipes**
- Individually wrapped
- Removes dirt and artifacts from your PSP photo exposure.





#### PROTECTIVE LEAD APRONS

Flow's lead aprons provide traditional protection for your patients in a variety of styles.

- · Considered the industry standard
- Manufactured from lead rubber with vinyl front
- Tear-resistant cling back
- · Easily cleaned with soap and water
- · Withstand years of heavy use



Thyroid collar 5inch x 27inch

22.00



Child 20inch x 21inch



Panoramic apron Adult 24inch x 27inch



Bib with attached collar 24inch x 26inch

# PREVENTATI AND HYGIENE

#### **FLEXIBRUSH**

Can your interdental brushes bend like this?



flexibility

With a 90° flexibility, Acers Flex Interdental Brushes allow you to access unreachable areas between your teeth.

- · High tensile strength
- · Easy grip handle
- 8 brushes per pack
- 24 packs/Carton
- · Price per pack

















#### PRODERMA DENT HAND MOISTURISER

#### Intensive Moisturising Hand Cream

- · Uniquely formulated for relief of dry hands due to latex gloves.
- Improves Seborrheic Dermatitis and Xerosis.
- · Revitalises skin.
- 125ml tube.





#### PRIME-DENT PRIME PASTE®



- · Prime-Paste is a blend of polishing and cleaning agents containing Sodium Fluoride
- This splatter-free formula rinses easily and completely
- · Medium grit 100g jar
- · Available in Raspberry, Bubblegum and Mint flavours.

#### PERFORATED **PROPHYLACTIC METAL STRIPS**

- 4mm wide. 100mm in length, 0.050mm thick
- Perforated prophylactic metal strips, recommended for tartar removal
- · 12pcs/packet.





#### IF-5 **BEAVERTAIL ULTRASONIC** SCALER **INSERT 25** AND 30 KHZ

- · This tip is best for removing hidden calculus
- · Compatible with most major brand magnetostrictive
- · CE and ISO certified and 510K registered
- Available in 25KHz and 30KHz.

# INFECTION CONTROL

#### NITRILE GLOVES

- · Exceptional comfort, fit and feel
- · Puncture and abrasion resistant
- · Chemical and oil resistant
- Puncture and abrasion resistant
- 100% nitrile rubber very skin friendly

OUR \$179.00 PER CARTON

- · Textured finger non-slip grip
- Exceptional comfort, fit and feel
- 240mm, 4.5g, 4.7mil at fingertip
- Dark blue
- Powder Free Nitrile.



#### LATEX DISPOSABLE GLOVES

- Exceptional comfort, fit and feel
- Puncture and abrasion resistant
- Powder free. Smooth.
- · 285MM, 10.5q
- Hand specific Left and Right hand
- 100 gloves/box 10 boxes per carton.









#### FFP2 **COMPLIANT**

· Classification FFP2 NR. EN 149:2001 + A1:2009 Respiratory **Protective Devices** - Filtering half masks to protect against particles -Requirements,



Testing, Marking. Non-valved. 5-layer mask. Nonsterile and disposable.

Certificate No: 2163-PPE-840.

OUR \$14.90



#### LEVEL 2 **SURGICAL FACE MASK**



- Fluid resistance of 120mmHg
- Elastic Earband (Earlooped)
- 3 Layers: Polypropylene spun-bond non-woven fabric (sanitary waterproof layer); Superfine polypropylene fibre melt-blown material layer; Polypropylene spunbond non-woven fabric (hygienic hygroscopic layer)
- Adjustable Nose Clip
- One Size Fits All Size: 175 x 95mm.



#### POLYPROPYLENE COATED IMPERVIOUS ISOLATION GOWN



- · Coated disposable isolation gowns provide effective fluid resistance from blood and body fluids whilst being breathable and comfortable. The disposable isolation gowns are designed for use with moderate to heavy fluids.
- Made from 20gm Polypropylene (PP) and 15g polyethylene (PE) the gowns are impervious on the front and sleeves.
- They have a velcro closure at the neck, full back with ties at each side and knitted cuffs.
- 1200mm x 1400mm x 40gsm.





#### SAFETY GLASSES

- Economical
- · Very lightweight, with thin streamline temples for better comfort
- · Lenses coated for excellent scratch resistance
- · Anti-fog coated, superior comfort and fit.









#### STERILIZED UNIVERSAL CONTACT LABELS

- · Latex free universal labels
- Can be used in Steam, EO, Gas and Hydrogen Peroxide sterilization
- 2 lines for Date and Load Number.





Code: SCL-010-K
OUR \$35.00

## STERILIZATION CONTACT LABEL GUN

 Lifetime warranty against defects or breakage under normal use conditions.



Code: SDG-001

OUR \$450.00

Code: IRC-001

OUR \$29.00

#### TYPE 6 STEAM INDICATOR

- Usage: Control of Steam sterilization processes.
- Emulation conditions: 3 minutes, 134 °C (FLASH cycles).
- Applicable Regulation: Designed under Quality Management Systems standards ISO 13485:2016/ NS-EN ISO 13485:2016.
- ISO 11140-1:2014
- Classification: Class 1, according to risk.
- Authorization: ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-4.
- Characteristics: Type 6 Emulator Indicator.
- One-side-laminated strips printed with indicator ink.
- Size: 70 mm x 23 mm.
- Initial color: yellow > Final color: black.





#### CHEMDYE® CD29



- Multi-variable indicator strips
- · Double card strips printed with indicator ink
- 200 x 15mm
- For steam sterilization processes
- · Turns from yellow to black
- Type 4 (ISO 11140-1).











#### STEAM INDICATOR STRIPS

- Class 1 Distinct change when processed
- 4 inch with date and operator writing space for permanent record.

4inch: SIS-250

OUR \$19.00

8inch: SIL-250
OUR \$29.00



## STEAM INDICATOR TAPE

- Colour changes when sterilised
- 55m rolls.

Blue linch: BT-036 | Yellow linch: ST-036 OUR PRICE \$19.00

#### STEAMPLUS STERILISATION INTEGRATOR



 Highest level of sterility assurance. Immediate Class 5 device for steam cycles. FDA approved as equivalent to Biological indicator.



#### TAPE DISPENSER

 Tape Dispenser can hold up to 2 rolls of 1 inch tape and has a weighted base.

Code: TD001
OUR \$69.00



#### CHEMDYE® HELIX - KIT

#### Helix PCD + Chemical Indicator Test Strips

- Pre-sterilization test kits for monitoring air-removal and steam penetration.
- Process Challenge Device for steam sterilization processes control
- For safe release of hollow loads with Type 2 chemical indicator
- Plastic capsule of great resistance and durability, to hold a chemical indicator designed for this system
- The capsule is connected to a high resistance plastic tube of 1.5m length.



Helix PCD: KH2X15

OUR \$249.00

PRICE



250 x replacements strips: KH2X15/PCD-A-18



- Chemical Indicator for Steam sterilization.
- Size: 125 x 125 x 16 mm.
- 134° C 3,5 min.
- Turns from Yellow to Dark Brown / Black.
- 100 % Lead Free.
- ISO 11140-1 Type 2 ISO 11140-4.

Code: BD125X/2

OUR \$169.00







#### PAPER TOWELS

#### 100% recycled paper content.



- · Recycled slimfold paper towels are a great economical choice for disposable paper hand towels. Made from 100% recycled tissue, the absorbent paper towels contain no dyes, inks or chlorine. FSC certified.
- N-fold sheet size of 230mm x 230mm
- Folded towel dimension 230mm x 76mm (1/3 Fold)
- 1 ply 200 sheet pack
- · Carton: 20 packets of 200 sheets (4000 sheets).

#### Reliable and durable cleaning.



- Slimfold paper towels are a standard, durable wiping solution for general use. Ideal for all cleaning requirements, quickly and easily cleaning all messes. FSC certified. Available in 1ply thickness.
- N-fold sheet size of 230mm x 230mm
- Folded towel dimension 230mm x 76mm (1/3 Fold)
- 1 ply 200 sheet pack
- · Carton: 20 packets of 200 sheets (4000 sheets).

OUR \$67.00

#### INCUBATOR DUAL FOR **BIOLOGICAL INDICATORS**

- · For biological indicators and hygiene monitoring systems 37± 2° or 60 ± 2°
- 26 Metal positions to incubate the BT20 indicators
- · Dimensions: 12cm high, 12.3cm diameter.





#### **BIOLOGICAL INDICATORS**

- · Self-contained Biological Indicator for Steam sterilization processes.Geobacillus Stearothermophilus (ATCC 7953) 10^5 or 10^6 spores per vial.
- · Monitoring vacuum-asissted and gravity air-displacement. Steam sterilization processes at 121-135°C. Readout = 24 hours
- Regulations: IRAM 37102-1, IRAM 37102-3, ISO 11138-1, ISO 11138-3.

#### RECORD NOTEBOOK

· Each notebook contains 17 pages.

#### A4 RECORD **NOTEBOOK**

· Each Sterilization Record Notebook comes with 100 pages.





Sterilizer



E sales@acers.co.nz

FREEPHONE **0800 422 377** 

w www.acers.co.nz